

NADL2 Polyclonal Antibody

Cat No: HR1AP12374

For research use only

Overview

Product Name	NADL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	87kDa
GeneID?Human?	NAALADL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>caution:Although related to the peptidase M28 family, it is probably inactive as an hydrolase due to lack of the conserved zinc-binding and active sites.,function:May be catalytically inactive.,miscellaneous:The gene maps to 3q26.31, a region associated with Cornelia de Lange syndrome. However, PubMed:15168106 failed to identify specific mutations in a panel of DNA samples from patients with Cornelia de Lange syndrome.,sequence caution:Contaminating sequence at the C-terminus.,similarity:Belongs to the peptidase M28 family. M28B subfamily.,tissue specificity:Expressed at higher level in kidney and placenta. In embryo, it is mainly confined to duodenal and stomach endoderm, mesonephros, metanephros and pancreas.,</p>