

## NALP2 Polyclonal Antibody

Cat No: HR1AP12369

For research use only

### Overview

Product Name	NALP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	116kDa
GeneID?Human?	NLRP2 NALP2 NBS1 PAN1 PYPAF2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NLR family pyrin domain containing 2(NLRP2) Homo sapiens NALP proteins, such as NALP2, are characterized by an N-terminal pyrin (MIM 608107) domain (PYD) and are involved in the activation of caspase-1 (CASP1; MIM 147678) by Toll-like receptors (see TLR4; MIM 603030). They may also be involved in protein complexes that activate proinflammatory caspases (Tschopp et al., 2003 [PubMed 12563287]).[supplied by OMIM, Mar 2008],