

CSPG4 Polyclonal Antibody

Cat No: HR1AP12365

For research use only

Overview

Product Name	CSPG4 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	255kDa
GeneID?Human?	CSPG4 MCSP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>chondroitin sulfate proteoglycan 4(CSPG4) Homo sapiens A human melanoma-associated chondroitin sulfate proteoglycan plays a role in stabilizing cell-substratum interactions during early events of melanoma cell spreading on endothelial basement membranes. CSPG4 represents an integral membrane chondroitin sulfate proteoglycan expressed by human malignant melanoma cells. [provided by RefSeq, Jul 2008],</p>