

## IGSF3 Polyclonal Antibody

Cat No: HR1AP12364

For research use only

### Overview

Product Name	IGSF3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	131kDa
GeneID?Human?	IGSF3 EWI3 KIAA0466
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	immunoglobulin superfamily member 3(IGSF3) Homo sapiens The protein encoded by this gene is an immunoglobulin-like membrane protein containing several V-type Ig-like domains. A mutation in this gene has been associated with bilateral nasolacrimal duct obstruction (LCDD). [provided by RefSeq, Jun 2016],