

NAAA Polyclonal Antibody

Cat No: HR1AP12361

For research use only

Overview

Product Name	NAAA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	NAAA AS AHL PLT
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	N-acylethanolamine acid amidase(NAAA) Homo sapiens This gene encodes an N-acylethanolamine-hydrolyzing enzyme which is highly similar to acid ceramidase. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],