

PAPP1 Polyclonal Antibody

Cat No: HR1AP12356

For research use only

Overview

Product Name	PAPP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	178kDa
GeneID?Human?	PAPPA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	pappalysin 1(PAPPA) Homo sapiens This gene encodes a secreted metalloproteinase which cleaves insulin-like growth factor binding proteins (IGFBPs). It is thought to be involved in local proliferative processes such as wound healing and bone remodeling. Low plasma level of this protein has been suggested as a biochemical marker for pregnancies with aneuploid fetuses. [provided by RefSeq, Jul 2008],