

DKK4 Polyclonal Antibody

Cat No: HR1AP12358

For research use only

Overview

Product Name	DKK4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	DKK4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	dickkopf WNT signaling pathway inhibitor 4(DKK4) Homo sapiens This gene encodes a protein that is a member of the dickkopf family. The secreted protein contains two cysteine rich regions and is involved in embryonic development through its interactions with the Wnt signaling pathway. Activity of this protein is modulated by binding to the Wnt co-receptor and the co-factor kremen 2. [provided by RefSeq, Jul 2008],