

PEAR1 Polyclonal Antibody

Cat No: HR1AP12355

For research use only

Overview

Product Name	PEAR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	114kDa
GeneID?Human?	PEAR1 MEGF12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	platelet endothelial aggregation receptor 1(PEAR1) Homo sapiens PEAR1 is a platelet receptor that signals upon the formation of platelet-platelet contacts independent of platelet activation and secondary to platelet aggregation (Nanda et al., 2005 [PubMed 15851471]).[supplied by OMIM, Mar 2008],