

PLXD1 Polyclonal Antibody

Cat No: HR1AP12347

For research use only

Overview

Product Name	PLXD1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	211kDa
GeneID?Human?	PLXND1 KIAA0620
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the plexin family.,similarity:Contains 1 Sema domain.,similarity:Contains 3 IPT/TIG domains.,tissue specificity:Detected at low levels in heart, placenta, lung, skeletal muscle, kidney, thymus and liver. Detected at very low levels in brain, colon, spleen, small intestine and peripheral blood leukocytes.,