

NDP Polyclonal Antibody

Cat No: HR1AP12348

For research use only

Overview

Product Name	NDP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	NDP EVR2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NDP, norrin cystine knot growth factor(NDP) Homo sapiens This gene encodes a secreted protein with a cystein-knot motif that activates the Wnt/beta-catenin pathway. The protein forms disulfide-linked oligomers in the extracellular matrix. Mutations in this gene result in Norrie disease and X-linked exudative vitreoretinopathy. [provided by RefSeq, Feb 2009],