

## PRS27 Polyclonal Antibody

Cat No: HR1AP12346

For research use only

### Overview

Product Name	PRS27 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	PRSS27 MPN UNQ1884/PRO4327
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>protease, serine 27(PRSS27) Homo sapiens This gene is located within a large protease gene cluster on chromosome 16. It belongs to the group-1 subfamily of serine proteases. The encoded protein is a secreted tryptic serine protease and is expressed mainly in the pancreas. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],</p>