

## ARC Polyclonal Antibody

Cat No: HR1AP12339

For research use only

### Overview

Product Name	ARC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	ARC KIAA0278
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Required for consolidation of synaptic plasticity as well as formation of long-term memory. Regulates endocytosis of AMPA receptors in response to synaptic activity. Required for homeostatic synaptic scaling of AMPA receptors.,similarity:Belongs to the ARC/ARG3.1 family.,subcellular location:Associated with the cell cortex of neuronal soma and dendrites. Enriched in postsynaptic density of dendritic spines. Associated with the sperm tail.,subunit:Interacts with SH3GL1/endophilin-2, SH3GL3/endophilin-3 and DNM2/DYN2.,</p>