

AMPO Polyclonal Antibody

Cat No: HR1AP12334

For research use only

Overview

Product Name	AMPO Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	90kDa
GeneID?Human?	AOPEP C9orf3 ONPEP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>chromosome 9 open reading frame 3(C9orf3) Homo sapiens This gene encodes a member of the M1 zinc aminopeptidase family. The encoded protein is a zinc-dependent metallopeptidase that catalyzes the removal of an amino acid from the amino terminus of a protein or peptide. This protein may play a role in the generation of angiotensin IV. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Oct 2010],</p>