

# CRTAM Polyclonal Antibody

Cat No: HR1AP12330

For research use only

## Overview

Product Name	CRTAM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	CRTAM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cytotoxic and regulatory T-cell molecule(CRTAM) Homo sapiens The CRTAM gene is upregulated in CD4 (see MIM 186940)-positive and CD8 (see CD8A; MIM 186910)-positive T cells and encodes a type I transmembrane protein with V and C1-like Ig domains (Yeh et al., 2008 [PubMed 18329370]).[supplied by OMIM, Feb 2009],