

TICN1 Polyclonal Antibody

Cat No: HR1AP12327

For research use only

Overview

Product Name	TICN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	SPOCK1 SPOCK TIC1 TICN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SPARC/osteonectin, cwcv and kazal like domains proteoglycan 1(SPOCK1) Homo sapiens This gene encodes the protein core of a seminal plasma proteoglycan containing chondroitin- and heparan-sulfate chains. The protein's function is unknown, although similarity to thyropin-type cysteine protease-inhibitors suggests its function may be related to protease inhibition. [provided by RefSeq, Jul 2008],