

ERMAP Polyclonal Antibody

Cat No: HR1AP12324

For research use only

Overview

Product Name	ERMAP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	ERMAP RD SC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	erythroblast membrane associated protein (Scianna blood group)(ERMAP) Homo sapiens The protein encoded by this gene is a cell surface transmembrane protein that may act as an erythroid cell receptor, possibly as a mediator of cell adhesion. Polymorphisms in this gene are responsible for the Scianna/Radin blood group system. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008],