

# LY86 Polyclonal Antibody

Cat No: HR1AP12320

For research use only

## Overview

Product Name	LY86 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	LY86 MD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:May cooperate with CD180 and TLR4 to mediate the innate immune response to bacterial lipopolysaccharide (LPS) and cytokine production. Important for efficient CD180 cell surface expression.,induction:In monocytes, down-regulated by the cell-wall fraction of Mycobacterium bovis (BCG-CWS).,subcellular location:Associated with CD180 at the cell surface.,subunit:Binds CD180.,tissue specificity:Highly expressed in B-cells, monocytes and tonsil.,</p>