

## S27A6 Polyclonal Antibody

Cat No: HR1AP12317

For research use only

### Overview

Product Name	S27A6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GeneID?Human?	SLC27A6 ACSVL2 FACVL2 FATP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 27 member 6(SLC27A6) Homo sapiens This gene encodes a member of the fatty acid transport protein family (FATP). FATPs are involved in the uptake of long-chain fatty acids and have unique expression patterns. Alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008],</p>