

HRK Polyclonal Antibody

Cat No: HR1AP12315

For research use only

Overview

Product Name	HRK Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	10kDa
GeneID?Human?	HRK BID3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>harakiri, BCL2 interacting protein(HRK) Homo sapiens This gene encodes a member of the BCL-2 protein family. Members of this family are involved in activating or inhibiting apoptosis. The encoded protein localizes to intracellular membranes. This protein promotes apoptosis by interacting with the apoptotic inhibitors BCL-2 and BCL-X(L) via its BH3 domain. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Oct 2012].</p>