

# NR6A1 Polyclonal Antibody

Cat No: HR1AP12310

For research use only

## Overview

Product Name	NR6A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	NR6A1 GCNF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nuclear receptor subfamily 6 group A member 1(NR6A1) Homo sapiens This gene encodes an orphan nuclear receptor which is a member of the nuclear hormone receptor family. Its expression pattern suggests that it may be involved in neurogenesis and germ cell development. The protein can homodimerize and bind DNA, but in vivo targets have not been identified. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Jun 2013],