

## ERR2 Polyclonal Antibody

Cat No: HR1AP12308

For research use only

### Overview

Product Name	ERR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	ESRRB ERRB2 ESRL2 NR3B2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	estrogen related receptor beta(ESRRB) Homo sapiens This gene encodes a protein with similarity to the estrogen receptor. Its function is unknown; however, a similar protein in mouse plays an essential role in placental development. [provided by RefSeq, Jul 2008].