

PAR15 Polyclonal Antibody

Cat No: HR1AP12309

For research use only

Overview

Product Name	PAR15 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	72kDa
GeneID?Human?	PARP15 BAL3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	poly(ADP-ribose) polymerase family member 15(PARP15) Homo sapiens PARP15 is a macrodomain-containing transcriptional repressor with poly(ADP-ribose) polymerase activity (Aguiar et al., 2005 [PubMed 16061477]).[supplied by OMIM, May 2008],