

# OTUD4 Polyclonal Antibody

Cat No: HR1AP12306

For research use only

## Overview

Product Name	OTUD4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	122kDa
GeneID?Human?	OTUD4 HIN1 KIAA1046
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	OTU deubiquitinase 4(OTUD4) Homo sapiens Alternatively spliced transcript variants have been found for this gene. The smaller protein isoform encoded by the shorter transcript variant is found only in HIV-1 infected cells. [provided by RefSeq, Jul 2010],