



CRIS3 Polyclonal Antibody

Cat No: HR1AP12299

For research use only

Overview

Product Name	CRIS3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	CRISP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the CRISP family.,subcellular location:In neutrophils, localized in specific granules.,tissue specificity:Salivary gland, pancreas and prostate > epididymis, ovary, thymus and colon.,