

CPXM1 Polyclonal Antibody

Cat No: HR1AP12289

For research use only

Overview

Product Name	CPXM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	CPXM1 CPX1 CPXM UNQ3015/PRO9782
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	carboxypeptidase X, M14 family member 1(CPXM1) Homo sapiens This gene likely encodes a member of the carboxypeptidase family of proteins. Cloning of a comparable locus in mouse indicates that the encoded protein contains a discoidin domain and a carboxypeptidase domain, but the protein appears to lack residues necessary for carboxypeptidase activity.[provided by RefSeq, May 2010],