

VSIG4 Polyclonal Antibody

Cat No: HR1AP12292

For research use only

Overview

Product Name	VSIG4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	VSIG4 CR1g Z39IG UNQ317/PRO362
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	V-set and immunoglobulin domain containing 4(VSIG4) Homo sapiens This gene encodes a v-set and immunoglobulin-domain containing protein that is structurally related to the B7 family of immune regulatory proteins. The encoded protein may be a negative regulator of T-cell responses. This protein is also a receptor for the complement component 3 fragments C3b and iC3b. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010],