

## BOC Polyclonal Antibody

Cat No: HR1AP12288

For research use only

### Overview

Product Name	BOC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	122kDa
GeneID?Human?	BOC UNQ604/PRO1190
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	BOC cell adhesion associated, oncogene regulated(BOC) Homo sapiens The protein encoded by this gene is a member of the immunoglobulin/fibronectin type III repeat family. It is a component of a cell-surface receptor complex that mediates cell-cell interactions between muscle precursor cells, and promotes myogenic differentiation. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2014],