

SG3A2 Polyclonal Antibody

Cat No: HR1AP12279

For research use only

Overview

Product Name	SG3A2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	10kDa
GeneID?Human?	SCGB3A2 PNSP1 UGRP1 UNQ566/PRO1128
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	secretoglobin family 3A member 2(SCGB3A2) Homo sapiens The protein encoded by this gene is a secreted lung surfactant protein and a downstream target of thyroid transcription factor. A single nucleotide polymorphism in the promoter of this gene results in susceptibility to asthma.[provided by RefSeq, Mar 2010],