

MCEM1 Polyclonal Antibody

Cat No: HR1AP12276

For research use only

Overview

Product Name	MCEM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	MCEMP1 C19orf59
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mast cell expressed membrane protein 1(MCEMP1) Homo sapiens This gene encodes a single-pass transmembrane protein. Based on its expression pattern, it is speculated to be involved in regulating mast cell differentiation or immune responses. [provided by RefSeq, Jul 2008],