

HES1 Polyclonal Antibody

Cat No: HR1AP12275

For research use only

Overview

Product Name	HES1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	HES1 BHLHB39 HL HRY
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	hes family bHLH transcription factor 1(HES1) Homo sapiens This protein belongs to the basic helix-loop-helix family of transcription factors. It is a transcriptional repressor of genes that require a bHLH protein for their transcription. The protein has a particular type of basic domain that contains a helix interrupting protein that binds to the N-box rather than the canonical E-box. [provided by RefSeq, Jul 2008],