

PAR14 Polyclonal Antibody

Cat No: HR1AP12273

For research use only

Overview

Product Name	PAR14 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	198kDa
GeneID?Human?	PARP14 BAL2 KIAA1268
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	poly(ADP-ribose) polymerase family member 14(PARP14) Homo sapiens This gene encodes a member of the poly(ADP-ribose) polymerase (PARP) protein family. The encoded anti-apoptotic protein may regulate aerobic glycolysis and promote survival of cancer cells. Increased expression of this gene has been reported in a variety of tumor types. [provided by RefSeq, Jul 2016],