

F19A2 Polyclonal Antibody

Cat No: HR1AP12271

For research use only

Overview

Product Name	F19A2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	FAM19A2 TAF2A2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>family with sequence similarity 19 member A2, C-C motif chemokine like(FAM19A2) Homo sapiens This gene is a member of the TAF2A family which is composed of five highly homologous genes that encode small secreted proteins. These proteins contain conserved cysteine residues at fixed positions, and are distantly related to MIP-1alpha, a member of the CC-chemokine family. The TAF2A proteins are predominantly expressed in specific regions of the brain, and are postulated to function as brain-specific chemokines or neurokines, that act as regulators of immune and nervous cells. [provided by RefSeq, Jul 2008],</p>