

NR1D2 Polyclonal Antibody

Cat No: HR1AP12268

For research use only

Overview

Product Name	NR1D2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	63kDa
GeneID?Human?	NR1D2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nuclear receptor subfamily 1 group D member 2(NR1D2) Homo sapiens This gene encodes a member of the nuclear hormone receptor family, specifically the NR1 subfamily of receptors. The encoded protein functions as a transcriptional repressor and may play a role in circadian rhythms and carbohydrate and lipid metabolism. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009],