

NR2E1 Polyclonal Antibody

Cat No: HR1AP12266

For research use only

Overview

Product Name	NR2E1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	NR2E1 TLX
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nuclear receptor subfamily 2 group E member 1(NR2E1) Homo sapiens The protein encoded by this gene is an orphan receptor involved in retinal development. The encoded protein also regulates adult neural stem cell proliferation and may be involved in control of aggressive behavior. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2015],