

NINJ1 Polyclonal Antibody

Cat No: HR1AP12257

For research use only

Overview

Product Name	NINJ1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	NINJ1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ninjurin 1(NINJ1) Homo sapiens The ninjurin protein is upregulated after nerve injury both in dorsal root ganglion neurons and in Schwann cells (Araki and Milbrandt, 1996 [PubMed 8780658]). It demonstrates properties of a homophilic adhesion molecule and promotes neurite outgrowth from primary cultured dorsal root ganglion neurons.[supplied by OMIM, Aug 2009],