

ATL1 Polyclonal Antibody

Cat No: HR1AP12254

For research use only

Overview

Product Name	ATL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	193kDa
GeneID?Human?	ADAMTSL1 ADAMTSR1 C9orf94 UNQ528/PRO1071
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ADAMTS like 1(ADAMTSL1) Homo sapiens This gene encodes a secreted protein and member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motif) family. This protein lacks the metalloproteinase and disintegrin-like domains, which are typical of the ADAMTS family, but contains other ADAMTS domains, including the thrombospondin type 1 motif. This protein may have important functions in the extracellular matrix. Alternative splicing results in multiple transcript variants encoding distinct proteins. [provided by RefSeq, Jul 2008].</p>