

GFRA3 Polyclonal Antibody

Cat No: HR1AP12255

For research use only

Overview

Product Name	GFRA3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	GFRA3 UNQ339/PRO538/PRO3664
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	GDNF family receptor alpha 3(GFRA3) Homo sapiens The protein encoded by this gene is a glycosylphosphatidylinositol(GPI)-linked cell surface receptor and a member of the GDNF receptor family. It forms a signaling receptor complex with RET tyrosine kinase receptor and binds the ligand, artemin. [provided by RefSeq, Jul 2008],