

ESAM Polyclonal Antibody

Cat No: HR1AP12253

For research use only

Overview

Product Name	ESAM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	ESAM UNQ220/PRO246
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Can mediate aggregation most likely through a homophilic molecular interaction.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,subunit:Interacts with BAIAP1.,tissue specificity:Highly expressed in endothelial cells.,