

C1QT5 Polyclonal Antibody

Cat No: HR1AP12248

For research use only

Overview

Product Name	C1QT5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	C1QTNF5 CTRP5 UNQ303/PRO344
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>C1q and tumor necrosis factor related protein 5(C1QTNF5) Homo sapiens This gene encodes a member of a family of proteins that function as components of basement membranes and may play a role in cell adhesion. Mutations in this gene have been associated with late-onset retinal degeneration. The protein may be encoded by either a bicistronic transcript including sequence from the upstream membrane frizzled-related protein gene (MFRP), or by a monocistronic transcript expressed from an internal promoter. [provided by RefSeq, Jun 2013],</p>