

DIAC Polyclonal Antibody

Cat No: HR1AP12247

For research use only

Overview

Product Name	DIAC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	CTBS CTB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	chitobiase(CTBS) Homo sapiens Chitobiase is a lysosomal glycosidase involved in degradation of asparagine-linked oligosaccharides on glycoproteins (Aronson and Kuranda, 1989 [PubMed 2531691]).[supplied by OMIM, Nov 2010],