

# LYVE1 Polyclonal Antibody

Cat No: HR1AP12242

For research use only

## Overview

Product Name	LYVE1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	LYVE1 CRSBP1 HAR XLKD1 UNQ230/PRO263
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	lymphatic vessel endothelial hyaluronan receptor 1(LYVE1) Homo sapiens This gene encodes a type I integral membrane glycoprotein. The encoded protein acts as a receptor and binds to both soluble and immobilized hyaluronan. This protein may function in lymphatic hyaluronan transport and have a role in tumor metastasis. [provided by RefSeq, Jul 2008],