

SIG14 Polyclonal Antibody

Cat No: HR1AP12244

For research use only

Overview

Product Name	SIG14 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	SIGLEC14 UNQ294/PRO333
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Putative adhesion molecule. Sialic acid-binding paired receptor which may activate associated receptors.,similarity:Belongs to the immunoglobulin superfamily. SIGLEC (sialic acid binding Ig-like lectin) family.,similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,similarity:Contains 2 Ig-like C2-type (immunoglobulin-like) domains.,subunit:Interacts with TYROBP.,tissue specificity:Mainly expressed in hematopoietic tissues including bone marrow, spleen and fetal liver. Also detected in lung and testis.,</p>