

ROBO4 Polyclonal Antibody

Cat No: HR1AP12240

For research use only

Overview

Product Name	ROBO4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	110kDa
GeneID?Human?	ROBO4 UNQ421/PRO3674
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Receptor for Slit proteins, at least for SLIT2, and seems to be involved in angiogenesis and vascular patterning. May mediate the inhibition of primary endothelial cell migration by Slit proteins.,similarity:Belongs to the immunoglobulin superfamily. ROBO family.,similarity:Contains 2 fibronectin type-III domains.,similarity:Contains 2 Ig-like C2-type (immunoglobulin-like) domains.,subunit:Interacts with SLIT2 and ENAH.,tissue specificity:Specifically expressed in endothelial cells. Expressed at sites of angiogenesis in different tumor types.,</p>