

TLN2 Polyclonal Antibody

Cat No: HR1AP12236

For research use only

Overview

Product Name	TLN2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	279kDa
GeneID?Human?	TLN2 KIAA0320
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>taln 2(TLN2) Homo sapiens This gene encodes a protein related to talin 1, a cytoskeletal protein that plays a significant role in the assembly of actin filaments and in spreading and migration of various cell types, including fibroblasts and osteoclasts. This protein has a different pattern of expression compared to talin 1 but, like talin 1, is thought to associate with unique transmembrane receptors to form novel linkages between extracellular matrices and the actin cytoskeleton. [provided by RefSeq, Jul 2008].</p>