

PYGO2 Polyclonal Antibody

Cat No: HR1AP12234

For research use only

Overview

Product Name	PYGO2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	PYGO2 PP7910
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Involved in signal transduction through the Wnt pathway.,sequence caution:Sequencing errors.,similarity:Contains 1 PHD-type zinc finger.,subunit:Binds to BCL9 via the PHD-type zinc finger motif, and thereby becomes part of the nuclear beta-catenin/TCF complex.,