

ANX13 Polyclonal Antibody

Cat No: HR1AP12231

For research use only

Overview

Product Name	ANX13 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	ANXA13 ANX13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>annexin A13(ANXA13) Homo sapiens This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The specific function of this gene has not yet been determined; however, it is associated with the plasma membrane of undifferentiated, proliferating endothelial cells and differentiated villus enterocytes. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008].</p>