

RGMA Polyclonal Antibody

Cat No: HR1AP12228

For research use only

Overview

Product Name	RGMA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	RGMA RGM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>repulsive guidance molecule family member a(RGMA) Homo sapiens This gene encodes a member of the repulsive guidance molecule family. The encoded protein is a glycosylphosphatidylinositol-anchored glycoprotein that functions as an axon guidance protein in the developing and adult central nervous system. This protein may also function as a tumor suppressor in some cancers. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009],</p>