

ANGL2 Polyclonal Antibody

Cat No: HR1AP12223

For research use only

Overview

Product Name	ANGL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	ANGPTL2 ARP2 UNQ170/PRO196
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>angiopoietin like 2(ANGPTL2) Homo sapiens Angiopoietins are members of the vascular endothelial growth factor family and the only known growth factors largely specific for vascular endothelium. Angiopoietin-1, angiopoietin-2, and angiopoietin-4 participate in the formation of blood vessels. ANGPTL2 protein is a secreted glycoprotein with homology to the angiopoietins and may exert a function on endothelial cells through autocrine or paracrine action. [provided by RefSeq, Jul 2008].</p>