

TCAM2 Polyclonal Antibody

Cat No: HR1AP12224

For research use only

Overview

Product Name	TCAM2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	TICAM2 TIRAP3 TIRP TRAM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	toll like receptor adaptor molecule 2(TICAM2) Homo sapiens TIRP is a Toll/interleukin-1 receptor (IL1R; MIM 147810) (TIR) domain-containing adaptor protein involved in Toll receptor signaling (see TLR4; MIM 603030).[supplied by OMIM, Apr 2004],