



FAM3A Polyclonal Antibody

Cat No: HR1AP12220

For research use only

Overview

Product Name	FAM3A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	FAM3A 2-19 2.19
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	family with sequence similarity 3 member A(FAM3A) Homo sapiens This gene encodes a cytokine-like protein. The expression of this gene may be regulated by peroxisome proliferator-activated receptor gamma, and the encoded protein may be involved in the regulation of glucose and lipid metabolism. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013],